




Part 1 Liter Vessel Inactinic Brown  
 Part Number 14143-12  
 Serial Number P36 000311 – P36 000360



# Test Certificate

		USP / Ph.Eur. Tolerances	
<b>Dimensions</b>	Vessel	Inside Ø a	102.9 to 105.1 mm
		Height b	168.2 mm
		With hemispherical bottom	160 – 210 mm
<b>Material</b>	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2024-03-21

21-03-2024 16:43:35

Moreshwar Pulekar  
 SOTAX Quality Management  
 C-64 | First Floor | Suchi Business Park | TTC Industrial area | MIDC Turbhe | Navi Mumbai 400705 | India